

## Product datasheet for **HP207842**

### NDUFB4 Human qPCR Primer Pair (NM\_004547)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4710
Forward Sequence:	GCTGAAACGAGAGTACCTGCTTC
Reverse Sequence:	TCCACACAGAGCTCCCATGAGT
ACCN:	<u><a href="#">NM_004547</a></u> , <u><a href="#">NM_004547.2</a></u> , <u><a href="#">NM_004547.3</a></u> , <u><a href="#">NM_004547.4</a></u> , <u><a href="#">NM_004547.5</a></u> , <u><a href="#">BC000855</a></u> , <u><a href="#">BC000855.1</a></u> , <u><a href="#">BC034579</a></u> , <u><a href="#">BC070285</a></u> , <u><a href="#">BC146885</a></u> , <u><a href="#">BC146887</a></u> , <u><a href="#">BC146927</a></u> , <u><a href="#">BC146932</a></u> , <u><a href="#">BF212059</a></u> , <u><a href="#">BF978138</a></u> , <u><a href="#">NM_004547.6</a></u>
UniProt ID:	<u><a href="#">O95168</a></u>
Synonyms:	B15; CI-B15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)