

Product datasheet for **HP207864**

PIK3C2G Human qPCR Primer Pair (NM_004570)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5288
Forward Sequence:	CTCCTTGCAGAGCATCCAGTT
Reverse Sequence:	GCACAGAATTGGAGAGCAGCTAG
ACCN:	<u>NM_004570</u> , <u>NM_004570.2</u> , <u>NM_004570.3</u> , <u>NM_004570.4</u> , <u>NM_004570.5</u> , <u>BC130277</u> , <u>BC144133</u> , <u>BC144134</u> , <u>BX092808</u>
UniProt ID:	<u>O75747</u>
Synonyms:	PI3K-C2-gamma; PI3K-C2GAMMA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)