

Product datasheet for **HP207866**

Phospholipase C beta 2 (PLCB2) Human qPCR Primer Pair (NM_004573)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5330
Forward Sequence:	CCTGGAAGTGACGGCTTATGAG
Reverse Sequence:	GCTCTGTGAAGGACGAGATGAC
ACCN:	NM_004573 , NM_004573.1 , NM_004573.2 , BC000939 , BC009009 , BC136467 , BT006905 , BU147302 , NM_004573.3
UniProt ID:	Q00722
Synonyms:	PLC-beta-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)