

Product datasheet for **HP207873**

RAB27A Human qPCR Primer Pair (NM_004580)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5873
Forward Sequence:	GAAGCCATAGCACTCGCAGAGA
Reverse Sequence:	CAGGACTTGTCCACACACCGTT
ACCN:	<u>NM_004580</u> , <u>NM_004580.1</u> , <u>NM_004580.2</u> , <u>NM_004580.3</u> , <u>NM_004580.4</u> , <u>BC132800</u> , <u>BC078139</u> , <u>BC107680</u> , <u>BC136423</u> , <u>BM014705</u> , <u>NM_004580.5</u>
UniProt ID:	<u>P51159</u>
Synonyms:	GS2; HsT18676; RAB27; RAM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)