

Product datasheet for **HP207892**

Syntaxin 1a (STX1A) Human qPCR Primer Pair (NM_004603)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6804
Forward Sequence:	TGGAGAACAGCATCCGTGAGCT
Reverse Sequence:	CCTCTCCACATAGTCTACCGCG
ACCN:	NM_004603 , NM_004603.1 , NM_004603.2 , NM_004603.3 , BC064644 , BC064644.1 , BC003011 , BC000444 , NM_004603.4
UniProt ID:	Q16623
Synonyms:	HPC-1; P35-1; STX1; SYN1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)