

Product datasheet for **HP207893**

Syntaxin 4 (STX4) Human qPCR Primer Pair (NM_004604)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6810
Forward Sequence:	CGGCAGACTATTGTCAAAGTGGG
Reverse Sequence:	ATCTCATCGCGCAGGTTCTGCA
ACCN:	NM_004604 , NM_004604.1 , NM_004604.2 , NM_004604.3 , NM_004604.4 , BC002436 , BC002436.2 , BI755734 , BQ775670 , BT007326 , NM_004604.5
UniProt ID:	Q12846
Synonyms:	p35-2; STX4A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)