

## Product datasheet for **HP207895**

### TAF1 Human qPCR Primer Pair (NM\_004606)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6872
Forward Sequence:	GGCTAAAGCTCTGCGCTGACTT
Reverse Sequence:	AGCACTGCTCTGGTGACACCAT
ACCN:	<u><a href="#">NM_004606</a></u> , <u><a href="#">NM_004606.1</a></u> , <u><a href="#">NM_004606.2</a></u> , <u><a href="#">NM_004606.3</a></u> , <u><a href="#">NM_004606.4</a></u> , <u><a href="#">BC172427</a></u> , <u><a href="#">BP284702</a></u> , <u><a href="#">BQ027170</a></u>
UniProt ID:	<u><a href="#">P21675</a></u>
Synonyms:	BA2R; CCG1; CCGS; DYT3; DYT3/TAF1; KAT4; MRXS33; N-TAF1; NSCL2; OF; P250; TAF(II)250; TAF2A; TAFII-250; TAFII250; XDP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)