

Product datasheet for **HP207911**

VIP Receptor 1 (VIPR1) Human qPCR Primer Pair (NM_004624)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7433
Forward Sequence:	GGACACCATCAACTCCTCACTG
Reverse Sequence:	CGCAGTTTCTGAAGCAGGATTCG
ACCN:	<u>NM_004624</u> , <u>NM_004624.1</u> , <u>NM_004624.2</u> , <u>NM_004624.3</u> , <u>BC064424</u> , <u>BC064424.1</u> , <u>BU681817</u> , <u>NM_004624.4</u>
UniProt ID:	<u>P32241</u>
Synonyms:	HVR1; II; PACAP-R-2; PACAP-R2; RDC1; V1RG; VAPC1; VIP-R-1; VIPR; VIRG; VPAC1; VPAC1R; VPCAP1R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)