

## Product datasheet for **HP207923**

### **UAP56 (DDX39B) Human qPCR Primer Pair (NM\_004640)**

#### **Product data:**

|                          |   |
|--------------------------|---|
| <b>Product Type:</b>     | qPCR Primer Pairs   |
| <b>Gene ID:</b>          | 7919  |
| <b>Forward Sequence:</b> | GATGTGCAGGATCGCTTTGAGG  |
| <b>Reverse Sequence:</b> | TCTCCTGAAGGACAGACGGTCA  |
| <b>ACCN:</b>             | <u><a href="#">NM_004640</a></u> , <u><a href="#">NM_004640.1</a></u> , <u><a href="#">NM_004640.2</a></u> , <u><a href="#">NM_004640.3</a></u> , <u><a href="#">NM_004640.4</a></u> , <u><a href="#">NM_004640.5</a></u> ,<br><u><a href="#">NM_004640.6</a></u> , <u><a href="#">BC013006</a></u> , <u><a href="#">BC013006.2</a></u> , <u><a href="#">BC000361</a></u> , <u><a href="#">BM993499</a></u> , <u><a href="#">BT009909</a></u> , <u><a href="#">BU628632</a></u> |
| <b>UniProt ID:</b>       | <u><a href="#">Q13838</a></u>   |
| <b>Synonyms:</b>         | BAT1; D6S81E; UAP56   |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| <b>Storage:</b>          | Store at -20°C.   |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)