

Product datasheet for **HP207957**

EIF1AY Human qPCR Primer Pair (NM_004681)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9086
Forward Sequence:	TGGAAGCATTGTGTTTTGATGGTG
Reverse Sequence:	CCTGATAGTCCCGTAGACCAAC
ACCN:	NM_004681 , NM_004681.1 , NM_004681.2 , NM_004681.3 , BC005248 , BC005248.1 , BF694967 , BG499107 , BP331805 , BT007201 , NM_004681.4
UniProt ID:	O14602
Synonyms:	eIF-4C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)