

Product datasheet for **HP207959**

SPARCL1 Human qPCR Primer Pair (NM_004684)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8404
Forward Sequence:	GTGAAGGCAACATGAGGGTGCA
Reverse Sequence:	GTTGGAGGACAAGTCACTGGATC
ACCN:	NM_004684 , NM_004684.1 , NM_004684.2 , NM_004684.3 , NM_004684.4 , NM_004684.5 , BC033721 , BC010383 , BC022484 , BX647713
UniProt ID:	Q14515
Synonyms:	MAST 9; MAST9; PIG33; SC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)