

Product datasheet for **HP207968**

MCT7 (SLC16A6) Human qPCR Primer Pair (NM_004694)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9120
Forward Sequence:	CGCTGTGTTTGCTTCGCACCA
Reverse Sequence:	TTTTCGGTGACGCTGGTCCTCT
ACCN:	NM_004694 , NM_004694.1 , NM_004694.2 , NM_004694.3 , NM_004694.4 , BC064832 , BC064832.1 , BC127749 , BC127750 , BC141962 , BC150547 , BC150548 , NM_004694.5
UniProt ID:	O15403
Synonyms:	MCT6; MCT7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)