

Product datasheet for **HP207971**

PRPF4 Human qPCR Primer Pair (NM_004697)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9128
Forward Sequence:	GGTAATGTGGCATCCTTCAGGAC
Reverse Sequence:	GATGCAGGATCTCCTCTTGAGC
ACCN:	NM_004697 , NM_004697.1 , NM_004697.2 , NM_004697.3 , NM_004697.4 , BC001588 , BC001588.2 , BC007424 , NM_004697.5
UniProt ID:	O43172
Synonyms:	HPRP4; HPRP4P; PRP4; Prp4p; RP70; SNRNP60
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)