

Product datasheet for **HP207980**

Apg12 (ATG12) Human qPCR Primer Pair (NM_004707)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9140
Forward Sequence:	GGAAGGACTTACGGATGTCTC
Reverse Sequence:	AGGAGTGTCTCCACAGCCTTT
ACCN:	<u>BC012266</u> , <u>NM_004707</u> , <u>NM_004707.1</u> , <u>NM_004707.2</u> , <u>NM_004707.3</u> , <u>BC012266.2</u> , <u>BC011033</u> , <u>BM697565</u> , <u>BP335591</u> , <u>BX327373</u>
UniProt ID:	<u>O94817</u>
Synonyms:	APG12; APG12L; FBR93; HAPG12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)