

## Product datasheet for **HP207983**

### Synaptogyrin 2 (SYNGR2) Human qPCR Primer Pair (NM\_004710)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9144
Forward Sequence:	GGTCATTGGTGACCTGCTCTTC
Reverse Sequence:	TGATGGCTGCCCTCACAGAGTC
ACCN:	<a href="#">NM_004710</a> , <a href="#">NM_004710.1</a> , <a href="#">NM_004710.2</a> , <a href="#">NM_004710.3</a> , <a href="#">NM_004710.4</a> , <a href="#">BC000407</a> , <a href="#">BC000407.2</a> , <a href="#">BC029755</a> , <a href="#">BC105992</a> , <a href="#">BF198109</a> , <a href="#">B1160076</a> , <a href="#">NM_004710.7</a>
UniProt ID:	<a href="#">O43760</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)