

Product datasheet for **HP207987**

DYRK1B Human qPCR Primer Pair (NM_004714)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9149
Forward Sequence:	GCCCGAAAACATCTTGCTGTGC
Reverse Sequence:	CACCTCAGGTGAGCGGTAGAAG
ACCN:	<u>BC025291</u> , <u>NM_004714</u> , <u>NM_004714.1</u> , <u>NM_004714.2</u> , <u>BC018751</u> , <u>BC018751.1</u> , <u>BC009688</u> , <u>BI964831</u> , <u>BQ446787</u> , <u>NM_004714.3</u>
UniProt ID:	<u>Q9Y463</u>
Synonyms:	AOMS3; MIRK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)