

Product datasheet for **HP207999**

DDX21 Human qPCR Primer Pair (NM_004728)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9188
Forward Sequence:	TGGACTCAGAGGGCAGCAGTTA
Reverse Sequence:	TGTCTCCATGCAAGGACTGAGC
ACCN:	<u>NM_004728</u> , <u>NM_004728.1</u> , <u>NM_004728.2</u> , <u>NM_004728.3</u> , <u>BC008071</u> , <u>BC008071.2</u> , <u>BC004182</u> , <u>BC009432</u> , <u>BC017228</u> , <u>BC104671</u> , <u>BU963561</u> , <u>BX648405</u> , <u>NM_004728.4</u>
UniProt ID:	<u>Q9NR30</u>
Synonyms:	GUA; GURDB; RH-II/GU; RH-II/GuA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)