

## Product datasheet for **HP208005**

### DCAMKL1 (DCLK1) Human qPCR Primer Pair (NM\_004734)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9201
Forward Sequence:	ACCGATGCCATCAAGCTGGACT
Reverse Sequence:	TCCTGGTAACGGAACCTTCTCCG
ACCN:	<u><a href="#">NM_004734</a></u> , <u><a href="#">NM_004734.1</a></u> , <u><a href="#">NM_004734.2</a></u> , <u><a href="#">NM_004734.3</a></u> , <u><a href="#">NM_004734.4</a></u> , <u><a href="#">BC152456</a></u> , <u><a href="#">BC160043</a></u> , <u><a href="#">BE672499</a></u> , <u><a href="#">BI824620</a></u> , <u><a href="#">BM679448</a></u> , <u><a href="#">BM682165</a></u> , <u><a href="#">BQ711413</a></u> , <u><a href="#">NM_004734.5</a></u>
UniProt ID:	<u><a href="#">O15075</a></u>
Synonyms:	CL1; CLICK1; DCAMKL1; DCDC3A; DCLK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)