

Product datasheet for **HP208010**

MTA2 Human qPCR Primer Pair (NM_004739)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9219
Forward Sequence:	TGAGATGGAGGAATGGTCAGCC
Reverse Sequence:	CTGGACTATGCTGGCAAGTGAC
ACCN:	NM_004739 , NM_004739.1 , NM_004739.2 , NM_004739.3 , BC053650 , BC053650.1 , BC023656 , BF590712 , BQ574987 , NM_004739.4
UniProt ID:	O94776
Synonyms:	MTA1L1; PID
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)