

Product datasheet for **HP208021**

DHRS3 Human qPCR Primer Pair (NM_004753)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9249
Forward Sequence:	GGGCACTGAGTGCCATTACTTC
Reverse Sequence:	CGGCATTGTTCCACGAGATGGT
ACCN:	NM_004753 , NM_004753.1 , NM_004753.2 , NM_004753.3 , NM_004753.4 , NM_004753.5 , NM_004753.6 , BC002730 , BC002730.2 , BF037805 , BX648476 , NM_004753.7
UniProt ID:	O75911
Synonyms:	DD83.1; RDH17; retSDR1; Rsd1; SDR1; SDR16C1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)