

## Product datasheet for **HP208040**

### MED1 Human qPCR Primer Pair (NM\_004774)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5469
Forward Sequence:	CTACAAAGGGCTGTCGGATGCA
Reverse Sequence:	GTGTCGGCTTGAATGGTTTCAGC
ACCN:	<u><a href="#">NM_004774</a></u> , <u><a href="#">NM_004774.1</a></u> , <u><a href="#">NM_004774.2</a></u> , <u><a href="#">BC006517</a></u> , <u><a href="#">BC060758</a></u> , <u><a href="#">BC131783</a></u> , <u><a href="#">BP357374</a></u> , <u><a href="#">NM_004774.4</a></u>
UniProt ID:	<u><a href="#">Q15648</a></u>
Synonyms:	CRSP1; CRSP200; DRIP205; DRIP230; PBP; PPARBP; PPARGBP; RB18A; TRAP220; TRIP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)