

## Product datasheet for **HP208042**

### **B4GALT5 Human qPCR Primer Pair (NM\_004776)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	9334
<b>Forward Sequence:</b>	GAAGATGACGACCTCTGGAACAG
<b>Reverse Sequence:</b>	GCCGTTCTTTTGACTTCCTCAGC
<b>ACCN:</b>	<u><a href="#">NM_004776</a></u> , <u><a href="#">NM_004776.1</a></u> , <u><a href="#">NM_004776.2</a></u> , <u><a href="#">NM_004776.3</a></u> , <u><a href="#">BC104987</a></u> , <u><a href="#">BC104987.1</a></u> , <u><a href="#">BC074821</a></u> , <u><a href="#">BC074873</a></u> , <u><a href="#">BC112265</a></u> , <u><a href="#">BG939833</a></u> , <u><a href="#">NM_004776.4</a></u>
<b>UniProt ID:</b>	<u><a href="#">O43286</a></u>
<b>Synonyms:</b>	B4Gal-T5; BETA4-GALT-IV; beta4Gal-T5; beta4GalT-V; gt-V
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)