

Product datasheet for **HP208061**

KIF3B Human qPCR Primer Pair (NM_004798)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9371
Forward Sequence:	ATCCTGGAGCAGAAACGACAGG
Reverse Sequence:	GTTCCAAGGTCTCCTCATCTCG
ACCN:	<u>NM_004798</u> , <u>NM_004798.1</u> , <u>NM_004798.2</u> , <u>NM_004798.3</u> , <u>BC136310</u> , <u>BC136311</u> , <u>BQ431652</u> , <u>BU621710</u> , <u>NM_004798.4</u>
UniProt ID:	<u>O15066</u>
Synonyms:	FLA8; HH0048; KLP-11; RP89
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)