

Product datasheet for **HP208072**

GRAP2 Human qPCR Primer Pair (NM_004810)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9402
Forward Sequence:	CCAAGAGGAGTGGTTTAAGGCG
Reverse Sequence:	GAGTAAGTTCTCTGCCTGGTGTC
ACCN:	<u>NM_004810</u> , <u>NM_004810.1</u> , <u>NM_004810.2</u> , <u>NM_004810.3</u> , <u>BC025692</u> , <u>BC025692.1</u> , <u>BC026002</u> , <u>NM_004810.4</u>
UniProt ID:	<u>O75791</u>
Synonyms:	GADS; GRAP-2; GRB2L; GRBLG; GrbX; Grf40; GRID; GRPL; Mona; P38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)