

## Product datasheet for **HP208073**

### Leupaxin (LPXN) Human qPCR Primer Pair (NM\_004811)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9404
Forward Sequence:	GAGATTGGCTCCAGTCCCTTCT
Reverse Sequence:	GGTCTGGTTCATTGCTGCAGC
ACCN:	<u><a href="#">NM_004811</a></u> , <u><a href="#">NM_004811.1</a></u> , <u><a href="#">NM_004811.2</a></u> , <u><a href="#">BC019035</a></u> , <u><a href="#">BC019035.2</a></u> , <u><a href="#">BC007490</a></u> , <u><a href="#">BC034230</a></u> , <u><a href="#">BP261314</a></u> , <u><a href="#">BU688891</a></u> , <u><a href="#">NM_004811.3</a></u>
UniProt ID:	<u><a href="#">O60711</a></u>
Synonyms:	LDPL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)