

Product datasheet for **HP208098**

HOMER2 Human qPCR Primer Pair (NM_004839)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9455
Forward Sequence:	CGTCACAGAAGTTTGGGCAGTG
Reverse Sequence:	CTTGGCAGCTTCTTTCACCTCC
ACCN:	NM_004839 , NM_004839.1 , NM_004839.2 , NM_004839.3 , BC012109 , BC012109.1 , BC044934 , BI666912 , NM_004839.4
UniProt ID:	Q9NSB8
Synonyms:	ACPD; CPD; DFNA68; HOMER-2; VESL-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)