

## Product datasheet for **HP208105**

### ATG5 Human qPCR Primer Pair (NM\_004849)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9474
Forward Sequence:	GCAGATGGACAGTTGCACACAC
Reverse Sequence:	GAGGTGTTCCAACATTGGCTCA
ACCN:	<a href="#">NM_004849</a> , <a href="#">NM_004849.1</a> , <a href="#">NM_004849.2</a> , <a href="#">NM_004849.3</a> , <a href="#">BC093011</a> , <a href="#">BC093011.1</a> , <a href="#">BC002699</a> , <a href="#">BM831472</a> , <a href="#">BP297166</a> , <a href="#">BX537904</a> , <a href="#">NM_004849.4</a>
UniProt ID:	<a href="#">Q9H1Y0</a>
Synonyms:	APG5; APG5-LIKE; APG5L; ASP; hAPG5; SCAR25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)