

Product datasheet for **HP208121**

SCAMP1 Human qPCR Primer Pair (NM_004866)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9522
Forward Sequence:	GAAACCAACAGAGGAACATCCAG
Reverse Sequence:	CCGACGATCTAATTCTGCGGCT
ACCN:	NM_004866 , NM_004866.1 , NM_004866.2 , NM_004866.3 , NM_004866.4 , NM_004866.5 , BC034048 , BC034048.1 , BC009787 , BC015065 , BC016355 , NM_004866.6
UniProt ID:	O15126
Synonyms:	SCAMP; SCAMP37
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)