

Product datasheet for **HP208128**

POLR1C Human qPCR Primer Pair (NM_004875)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9533
Forward Sequence:	GTTTCGTCTCCAGGTCAGATGC
Reverse Sequence:	CCTCTGGAAAGAGATCAGCCTG
ACCN:	<u>BC008863</u> , <u>NM_004875</u> , <u>NM_004875.1</u> , <u>NM_004875.2</u> , <u>BC008863.2</u> , <u>BC008118</u> , <u>BM828768</u> , <u>BQ575535</u>
UniProt ID:	<u>O15160</u>
Synonyms:	RP3-337H4.4; RPA5; RPA39; RPA40; RPAC1; TCS3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)