

Product datasheet for **HP208142**

SEC22L1 (SEC22B) Human qPCR Primer Pair (NM_004892)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9554
Forward Sequence:	GATGTGCAGAGGATCATGGTGG
Reverse Sequence:	CTTCGCATCCTGGCGGTATTTTC
ACCN:	<u>NM_004892</u> , <u>NM_004892.1</u> , <u>NM_004892.2</u> , <u>NM_004892.3</u> , <u>NM_004892.4</u> , <u>NM_004892.5</u> , <u>BC001364</u> , <u>BC001364.2</u> , <u>BU741516</u>
UniProt ID:	<u>O75396</u>
Synonyms:	ERS-24; SEC22L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)