

Product datasheet for **HP208189**

ER81 (ETV1) Human qPCR Primer Pair (NM_004956)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2115
Forward Sequence:	GCAAGAAGGCTTCCTGGCTCAT
Reverse Sequence:	CCTCCCGATACATTCCTGGCT
ACCN:	NM_004956 , NM_004956.1 , NM_004956.2 , NM_004956.3 , NM_004956.4 , BC098403 , BC098403.1 , BC042950 , BC045776 , BC106762 , BC106763 , BF346325 , BU737178 , NM_004956.5
UniProt ID:	P50549
Synonyms:	ER81
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)