

Product datasheet for **HP208213**

LIMS1 Human qPCR Primer Pair (NM_004987)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3987
Forward Sequence:	TGGTGCGTGAAGTCTTTGCCT
Reverse Sequence:	CATAGCACTTCTTACAGACTGGC
ACCN:	NM_004987 , NM_004987.1 , NM_004987.2 , NM_004987.3 , NM_004987.4 , NM_004987.5 , BC005341 , BC005341.1 , BC015843 , BC022365 , BG942518 , NM_004987.6
UniProt ID:	P48059
Synonyms:	PINCH; PINCH-1; PINCH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)