

Product datasheet for **HP208216**

EVI1 (MECOM) Human qPCR Primer Pair (NM_004991)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2122
Forward Sequence:	CCTTATGTGGGAGAGCAGAGGT
Reverse Sequence:	GAAGGCTATTCTACGTCTGAGC
ACCN:	<u>NM_004991</u> , <u>NM_004991.3</u> , <u>BC069498</u> , <u>BC143952</u> , <u>BC130520</u> , <u>BC012450</u> , <u>BC031019</u> , <u>BC143951</u> , <u>BX484790</u> , <u>BX640908</u> , <u>BX647613</u> , <u>NM_004991.4</u>
UniProt ID:	<u>Q13465</u>
Synonyms:	AML1-EVI-1; EVI1; KMT8E; MDS1; MDS1-EVI1; PRDM3; RUSAT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)