

Product datasheet for **HP208218**

MMP14 Human qPCR Primer Pair (NM_004995)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4323
Forward Sequence:	CCTTGACTGTCAGGAATGAGG
Reverse Sequence:	TTCTCCGTGCCATCCACTGGT
ACCN:	<u>NM_004995</u> , <u>NM_004995.1</u> , <u>NM_004995.2</u> , <u>NM_004995.3</u> , <u>BC064803</u> , <u>BC064803.1</u> , <u>BC039581</u> , <u>BQ009340</u> , <u>NM_004995.4</u>
UniProt ID:	<u>P50281</u>
Synonyms:	MMP-14; MMP-X1; MT-MMP; MT-MMP 1; MT1-MMP; MT1MMP; MTMMP1; WNCHRS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)