

Product datasheet for **HP208239**

Profilin 1 (PFN1) Human qPCR Primer Pair (NM_005022)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5216
Forward Sequence:	CATCGTGGGCTACAAGGACTCG
Reverse Sequence:	CCAAGTGTGAGCCATTCACGT
ACCN:	NM_005022 , NM_005022.1 , NM_005022.2 , NM_005022.3 , BC015164 , BC015164.1 , BC002475 , BC006768 , BC013439 , BC057828 , BM850671 , BT007001 , NM_005022.4
UniProt ID:	P07737
Synonyms:	ALS18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)