

Product datasheet for **HP208244**

PIP5K2 alpha (PIP4K2A) Human qPCR Primer Pair (NM_005028)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5305
Forward Sequence:	GCCAGATGACTTCAAAGCCTATTC
Reverse Sequence:	CCGCAGGTTACGGAAGACCATC
ACCN:	<u>NM_005028</u> , <u>NM_005028.1</u> , <u>NM_005028.2</u> , <u>NM_005028.3</u> , <u>NM_005028.4</u> , <u>BC018034</u> , <u>BC018034.1</u> , <u>NM_005028.5</u>
UniProt ID:	<u>P48426</u>
Synonyms:	PI5P4KA; PIP5K2A; PIP5KII-alpha; PIP5KIIA; PIPK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)