

Product datasheet for **HP208252**

Cyclophilin 40 (PPID) Human qPCR Primer Pair (NM_005038)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5481
Forward Sequence:	CAGAATGGGACAGGTGGAGAAAG
Reverse Sequence:	CTGAGAACCGTTTGTGTTGCGG
ACCN:	<u>NM_005038</u> , <u>NM_005038.1</u> , <u>NM_005038.2</u> , <u>BC030707</u> , <u>BC030707.2</u> , <u>BC014893</u> , <u>BG699755</u> , <u>BU726810</u> , <u>NM_005038.3</u>
UniProt ID:	<u>Q08752</u>
Synonyms:	CYP-40; CYPD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)