

## Product datasheet for **HP208263**

### QARS1 Human qPCR Primer Pair (NM\_005051)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5859
Forward Sequence:	CACACTACGGATGAAGCTGGTG
Reverse Sequence:	GATGGAGTCACAGAGGCAGTGT
ACCN:	<u><a href="#">NM_005051</a></u> , <u><a href="#">NM_005051.1</a></u> , <u><a href="#">NM_005051.2</a></u> , <u><a href="#">BC000394</a></u> , <u><a href="#">BC000394.2</a></u> , <u><a href="#">BC001567</a></u> , <u><a href="#">BC001772</a></u> , <u><a href="#">BC016634</a></u> , <u><a href="#">BC029739</a></u> , <u><a href="#">BM560535</a></u> , <u><a href="#">NM_005051.3</a></u>
UniProt ID:	<u><a href="#">P47897</a></u>
Synonyms:	GLNRS; MSCCA; PRO2195; QARS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)