

Product datasheet for **HP208280**

SLC12A4 Human qPCR Primer Pair (NM_005072)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6560
Forward Sequence:	CTGTAACAGGCATCATGGCTGG
Reverse Sequence:	CACACTGCTGAAGTACACGAGG
ACCN:	NM_005072 , NM_005072.1 , NM_005072.2 , NM_005072.3 , NM_005072.4 , BC021193 , BC021193.2 , BC035480
UniProt ID:	Q9UP95
Synonyms:	CTC-479C5.17; hKCC1; KCC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)