

Product datasheet for **HP208303**

ISG15 Human qPCR Primer Pair (NM_005101)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9636
Forward Sequence:	CTCTGAGCATCCTGGTGAGGAA
Reverse Sequence:	AAGGTCAGCCAGAACAGGTCGT
ACCN:	<u>BC009507</u> , <u>NM_005101</u> , <u>NM_005101.1</u> , <u>NM_005101.2</u> , <u>NM_005101.3</u> , <u>BC009507.2</u> , <u>BM712238</u> , <u>BT007297</u> , <u>BU161226</u> , <u>NM_005101.4</u>
UniProt ID:	<u>P05161</u>
Synonyms:	G1P2; hUCRP; IFI15; IMD38; IP17; UCRP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)