

Product datasheet for **HP208307**

ENDOGL1 (EXOG) Human qPCR Primer Pair (NM_005107)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9941
Forward Sequence:	CCTCAGACTAGAGGCGATGGAA
Reverse Sequence:	CGCCAGTGGTTCGGTAGATACT
ACCN:	NM_005107 , NM_005107.1 , NM_005107.2 , NM_005107.3 , BC015820 , BC065368 , BC104212 , BC104213 , NM_005107.4
UniProt ID:	Q9Y2C4
Synonyms:	ENDOGL1; ENDOGL2; ENGL; ENGL-a; ENGL-b; ENGLA; ENGLB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)