

Product datasheet for **HP208310**

GFPT2 Human qPCR Primer Pair (NM_005110)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9945
Forward Sequence:	GCTCATCGTGATTGGCTGTGGA
Reverse Sequence:	CAACCATCACAGGAAGCTCAGTC
ACCN:	<u>NM_005110</u> , <u>NM_005110.1</u> , <u>NM_005110.2</u> , <u>NM_005110.3</u> , <u>BC000012</u> , <u>BC000012.2</u> , <u>BI753491</u> , <u>BT009818</u> , <u>NM_005110.4</u>
UniProt ID:	<u>O94808</u>
Synonyms:	GFAT; GFAT 2; GFAT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)