

## Product datasheet for **HP208318**

### TL1A (TNFSF15) Human qPCR Primer Pair (NM\_005118)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9966
Forward Sequence:	CACCACATACCTGCTTGTCAGC
Reverse Sequence:	TCTCCGTCTGCTCTAAGAGGTG
ACCN:	<u><a href="#">NM_005118</a></u> , <u><a href="#">NM_005118.1</a></u> , <u><a href="#">NM_005118.2</a></u> , <u><a href="#">NM_005118.3</a></u> , <u><a href="#">BC104463</a></u> , <u><a href="#">BC069435</a></u> , <u><a href="#">BC074940</a></u> , <u><a href="#">BC074941</a></u> , <u><a href="#">BC104462</a></u> , <u><a href="#">NM_005118.4</a></u>
UniProt ID:	<u><a href="#">O95150</a></u>
Synonyms:	TL1; TL1A; TNLG1B; VEGI; VEGI192A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)