

Product datasheet for **HP208330**

DGCR2 Human qPCR Primer Pair (NM_005137)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9993
Forward Sequence:	GTGCTGCAAGTTCATGTGTCTGG
Reverse Sequence:	CAGCAGTGACAGGATGAGGAAG
ACCN:	NM_005137 , NM_005137.2 , BC040500 , BC040500.1 , BC020339 , BC032430 , BC035425 , NM_005137.3
UniProt ID:	P98153
Synonyms:	DGS-C; IDD; LAN; SEZ-12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)