

Product datasheet for **HP208331**

SCO2 Human qPCR Primer Pair (NM_005138)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9997
Forward Sequence:	GACCACTCCATTGCCATCTACC
Reverse Sequence:	CTCAAGACAGGACACTGCGGAA
ACCN:	<u>NM_005138</u> , <u>NM_005138.1</u> , <u>NM_005138.2</u> , <u>BC102024</u> , <u>BC102025</u> , <u>BG392981</u> , <u>BM787182</u> , <u>BQ878572</u> , <u>NM_005138.3</u>
UniProt ID:	<u>O43819</u>
Synonyms:	CEMCOX1; ECGF1; Gliostatin; MC4DN2; MYP6; PD-ECGF; SCO1L; TdRPase; TP; TYMP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)