

Product datasheet for **HP208336**

Haptoglobin (HP) Human qPCR Primer Pair (NM_005143)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3240
Forward Sequence:	TGGCTATGTGGAGCACTCGGTT
Reverse Sequence:	CAGGAAGTTTATCTCCAACAGCC
ACCN:	<u>BC058031</u> , <u>NM_005143</u> , <u>NM_005143.1</u> , <u>NM_005143.2</u> , <u>NM_005143.3</u> , <u>NM_005143.4</u> , <u>BC058031.1</u> , <u>BC107587</u> , <u>BC017862</u> , <u>BC070299</u> , <u>BC121124</u> , <u>BC121125</u> , <u>BJ994335</u> , <u>NM_005143.5</u>
UniProt ID:	<u>P00738</u>
Synonyms:	BP; HP2ALPHA2; HPA1S
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)