

Product datasheet for **HP208359**

ATOH1 Human qPCR Primer Pair (NM_005172)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	474
Forward Sequence:	CCTTCCAGCAAACAGGTGAATGG
Reverse Sequence:	GAACGACGGGATAACATTGCGC
ACCN:	NM_005172 , NM_005172.1 , BC069594 , BC069594.1 , BC069145 , BC069578 , BC069604 , BC113623 , BC113625 , NM_005172.2
UniProt ID:	Q92858
Synonyms:	ATH1; bHLHa14; HATH1; MATH-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)