

## Product datasheet for **HP208366**

### CALM3 Human qPCR Primer Pair (NM\_005184)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	808
Forward Sequence:	GAGAGCGTTCCGTGTCTTTGA
Reverse Sequence:	ACCTCCTCATCGGTCAGCTTCT
ACCN:	<a href="#">NM_005184</a> , <a href="#">NM_005184.1</a> , <a href="#">NM_005184.2</a> , <a href="#">NM_005184.3</a> , <a href="#">BC005137</a> , <a href="#">BC006182</a> , <a href="#">BG753029</a> , <a href="#">BQ428612</a> , <a href="#">BT006855</a> , <a href="#">NM_005184.4</a>
UniProt ID:	<a href="#">P62158</a>
Synonyms:	CALM; CaM; CAM1; CAM2; CAMB; CaMIII; CPVT6; HEL-S-72; LQT16; PHKD; PHKD3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)