

Product datasheet for **HP208379**

CCR8 Human qPCR Primer Pair (NM_005201)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1237
Forward Sequence:	TGGCTGTTGTCCATGCCGTGTA
Reverse Sequence:	TGGGATGGTAGCCATAATGGCG
ACCN:	NM_005201 , NM_005201.1 , NM_005201.2 , NM_005201.3 , BC069067 , BC074743 , BC107152 , BC107159 , NM_005201.4
UniProt ID:	P51685
Synonyms:	CC-CKR-8; CCR-8; CDw198; CKRL1; CMKBR8; CMKBRL2; CY6; GPRCY6; TER1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)