

Product datasheet for **HP208393**

DDOST Human qPCR Primer Pair (NM_005216)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 1650 |
| Forward Sequence: | TACGCTCATCGTGGCTGACACT |
| Reverse Sequence: | CCAGCACCAAAGGTTATCAGG |
| ACCN: | <u>NM_005216</u> , <u>NM_005216.1</u> , <u>NM_005216.2</u> , <u>NM_005216.3</u> , <u>NM_005216.4</u> , <u>BC002594</u> , <u>BC002594.2</u> , <u>BP233191</u> |
| UniProt ID: | <u>P39656</u> |
| Synonyms: | AGER1; CDG1R; GATD6; OKSWcl45; OST; OST48; WBP1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)